






WESTAR substrates for WB

WESTAR substrates are Cyanagen's patented product line for Western Blotting, which offer several advantages over other detection reagents.

- A broad range of sensitivities
- Long light emission
- Immediate technical support

	Primary benefits	Detection level	Recommended Ab dilution								
WESTAR ECL-Sun											
 <p>1° Ab 1/1,000 - 2° Ab 1/15,000</p>	Direct alternative to the Pierce ECL with more stable working solution (5 days at 4°C)	~ 20 picograms	<table border="0"> <tr> <td>1°</td> <td>1/100</td> </tr> <tr> <td></td> <td>1/5,000</td> </tr> <tr> <td>2°</td> <td>1/1,000</td> </tr> <tr> <td></td> <td>1/15,000</td> </tr> </table>	1°	1/100		1/5,000	2°	1/1,000		1/15,000
1°	1/100										
	1/5,000										
2°	1/1,000										
	1/15,000										
WESTAR Nova 2011											
 <p>1° Ab 1/2,000 - 2° Ab 1/30,000</p>	Provides excellent performance for routine experiments	~ 1 picogram	<table border="0"> <tr> <td>1°</td> <td>1/500</td> </tr> <tr> <td></td> <td>1/5,000</td> </tr> <tr> <td>2°</td> <td>1/20,000</td> </tr> <tr> <td></td> <td>1/100,000</td> </tr> </table>	1°	1/500		1/5,000	2°	1/20,000		1/100,000
1°	1/500										
	1/5,000										
2°	1/20,000										
	1/100,000										
WESTAR EtaC											
 <p>1° Ab 1/4,000 - 2° Ab 1/40,000</p>	High signal to noise ratio for medium-abundance proteins	~ 500 femtograms	<table border="0"> <tr> <td>1°</td> <td>1/1,000</td> </tr> <tr> <td></td> <td>1/15,000</td> </tr> <tr> <td>2°</td> <td>1/25,000</td> </tr> <tr> <td></td> <td>1/150,000</td> </tr> </table>	1°	1/1,000		1/15,000	2°	1/25,000		1/150,000
1°	1/1,000										
	1/15,000										
2°	1/25,000										
	1/150,000										
WESTAR EtaC Ultra											
 <p>1° Ab 1/20,000 - 2° Ab 1/100,000</p>	Very sensitive with extended signal duration for long exposure	~ 150 femtograms	<table border="0"> <tr> <td>1°</td> <td>1/5,000</td> </tr> <tr> <td></td> <td>1/50,000</td> </tr> <tr> <td>2°</td> <td>1/50,000</td> </tr> <tr> <td></td> <td>1/250,000</td> </tr> </table>	1°	1/5,000		1/50,000	2°	1/50,000		1/250,000
1°	1/5,000										
	1/50,000										
2°	1/50,000										
	1/250,000										
WESTAR Supernova											
 <p>1° Ab 1/30,000 - 2° Ab 1/200,000</p>	Most sensitive substrate for high-end blotting	~ 50 femtograms	<table border="0"> <tr> <td>1°</td> <td>1/5,000</td> </tr> <tr> <td></td> <td>1/100,000</td> </tr> <tr> <td>2°</td> <td>1/100,000</td> </tr> <tr> <td></td> <td>1/500,000</td> </tr> </table>	1°	1/5,000		1/100,000	2°	1/100,000		1/500,000
1°	1/5,000										
	1/100,000										
2°	1/100,000										
	1/500,000										

Detection of serial dilutions of human transferrin (from 5.0 ng to 0.5 ng) performed using Cyanagen WESTAR chemiluminescent substrates. The blots were developed using anti-transferrin antibody and goat anti-rabbit HRP at the indicated dilutions.

Your local distributor:



THE INNOVATOR FOR SCIENCE

WWW.ISOGEN-LIFESCIENCE.COM

Isogen Life Science B.V.
Veldzigt 2a
3454 PW De Meern
The Netherlands
tel. +31 30 688 0771

For service questions, please contact
service@isogen-lifescience.com

For all other questions, please contact
support@isogen-lifescience.com