

# CEInfinite

Accelerate Your Biologics Development and  
Protein Characterization




The Next Generation of iCIEF

## The CEInfinite Instrument


### Imaged Capillary Isoelectric Focusing (iCIEF)

The CEInfinite system is based on proprietary whole column imaging detection (WCID) technology. A pH gradient is created using an electric field and carrier ampholytes, allowing proteins to separate along the capillary according to their isoelectric points (pIs). Highly efficient separation is achieved when proteins are focused into narrow zones in the pH gradient. An advanced scientific CMOS (sCMOS) image sensor and a deep UV LED light source is used to monitor real time separation in the column.




### High Separation Efficiency

In comparison to traditional CIEF and gel electrophoresis, proteins separated with the CEInfinite instrument do not need mobilization and post-run staining steps which greatly shortens sample analysis time. Each run takes approximately 10 minutes, streamlining method development from months to hours and separation resolution can be as low as 0.01 pH.



## Protein Charge Heterogeneity Analysis

Protein charge heterogeneity profiling and characterization is essential to biologics development. This assay is performed quickly and easily, with little method development required, with the CEInfinite system. Results are generated with great reproducibility not only from within the same sample vial, but also across several sample vials, ensuring precision during both analysis and fractionation.



	pI		Peak Area	
	Average	% RSD	Average	% RSD
A1	9.18	0.06	32.67	1.63
mAb	9.23	0.00	59.62	0.82
B1	9.31	0.07	6.83	2.56
B2	9.41	0.05	0.92	13.42


iCIEF of NIST mAb taken from six different vials and the distribution of charge profiles

## Total Solutions


Support for your CEInfinite instrument is available and includes novel applications such as high-resolution iCIEF fractionation, iCIEF coupling to MS, protein interaction studies, and Taylor dispersion rapid measurement of protein hydrodynamic radius.

## Isolate Products for Protein Characterization (Patent Pending)

The high resolution of the CEInfinite iCIEF system translates to the utmost purity of the final protein product. For example, proteins with a difference in pI of 0.2 or greater can be fractionated with greater than 80% purity. With the new CEInfinite 320  $\mu\text{m}$  cartridge, one 45-minute run can yield up to 40  $\mu\text{g}$  of purified protein.




iCIEF separation of mAb charge variants



Proprietary cartridges are designed to prevent peak remixing

## iCIEF Separation of Antibody-Drug Conjugates

The CEInfinite system also provides high-resolution separation of antibody-drug conjugates (ADCs).



## Protein Molecular Characterization by iCIEF-MS

CEInfinite cartridges for iCIEF-MS preserve the great separation resolution of iCIEF. The use of a proprietary solvent during separation eliminates the need for urea and makes the protein fractions produced ready for MS with the interfacing technology available to connect the iCIEF-MS cartridge to an ESI source for MS.

## CEInfinite Base, Analytical, and Preparative Systems

Three CEInfinite systems are currently available, each with high-performance iCIEF capabilities.

	Base	Analytical	Preparative
iCIEF Separation	✓	✓	✓
Whole Column Image Detection	✓	✓	✓
Autosampler	✗	✓	✓
Pump	✗	✗	✓
Fractionation Capabilities	✗	✗	✓
iCIEF Coupled to Mass Spectrometry	✗	✗	✓

## iCIEF Complete Consumable Product Line

### Cartridges



CEInfinite cartridges are the world's first WCID imaging cartridges for direct coupling to a mass spectrometer. In addition to fluorocarbon (FC) and polyacrylamide (PA), proprietary acrylamide derivatives (AD) provide a highly reliable molecular layer coating suitable for separations between pH 2-11. These cartridges deliver high sensitivity and performance, and the easy flow along the capillary makes for simpler, faster sample injection with less clogging. Cartridges are available with various inner diameters for increased sensitivity and capacity.

### AESlytes



AESlytes are synthesized strictly to ensure a high performance CIEF separation. Their low interaction with proteins increases resolving power and AESlytes do not strongly absorb UV, eliminating baseline noise. Finally, their linear pH gradient makes for a remarkably accurate pI calibration.

### pI Markers



Reference standard isoelectric point (pI) markers are essential for calibrating your iCIEF. The CEInfinite system brings you the most comprehensive, evenly spaced, and thoroughly calibrated pI markers in the world. All CEInfinite pI markers come in a ready-to-use 100X concentration and have long shelf lives, reducing waste and saving you restocking time.

### Electrolytes & Polymer Solutions



Manufactured under strict quality control guidelines, CEInfinite consumables will save you time as well and eliminate variability for a more reproducible result. All CEInfinite test kits, electrolytes, and methylcellulose solutions are fully tested and certified.

## Product Specifications

### iCIEF Instrument

<b>Detection Mode</b>	Whole column, sCMOS imaging technology
<b>Height x Width x Depth</b>	54 cm x 33 cm x 30 cm
<b>Weight</b>	30 pounds (14 kg)
<b>Detection Dynamic Range</b>	250 (0.004 – 1.0 AU, 280 nm)
<b>High Voltage Range</b>	0 – 3000 Volts (Continuously Adjustable)
<b>Sample Throughput</b>	Up to 12 injections per hour
<b>Working Temperature</b>	15 - 35°C
<b>Humidity</b>	20 - 80% RH
<b>Electrical Requirement</b>	100/240 VAC, 50 - 60 Hz
<b>Exposure Time</b>	0.02 – 99.9 ms
<b>ADC Maximum</b>	16386 AU
<b>Operation Mode</b>	Manual or Automatic
<b>Detector Noise</b>	Less than 0.001 AU, 280 nm
<b>Separation pH Range</b>	2.1 – 11.1

### Autosampler

<b>Model</b>	840
<b>Sample Capacity</b>	84+3 vial tray, 96 well plate
<b>Sample Tray Temperature</b>	4 - 40°C
<b>Typical Sample Volume</b>	15 µL
<b>Electrical Requirement</b>	95 – 240 VAC, 50 – 60 Hz
<b>Height x Width x Depth</b>	36 cm x 30 cm x 57.5 cm
<b>Weight</b>	46 pounds (21 kg)
<b>Working Temperature</b>	10 - 40°C
<b>Humidity</b>	20 – 80% RH

### Syringe Pump

<b>Syringe Size</b>	250 µL
<b>Electrical Requirement</b>	100 – 240 VAC, 50 – 60 Hz
<b>Height x Width x Depth</b>	12.14 cm x 10.8 cm x 24.1 cm
<b>Weight</b>	2.72 pounds (1.23kg)
<b>Working Temperature</b>	5 – 40°C
<b>Humidity</b>	20 – 80% RH

### Software

CEInsight software fully controls the CEInfinite system with a user-friendly GUI. The CFR control and quantitation software provides reliable operation control, highly efficient data processing, and FDA Title 21 CFR Part 11 compliance.



## **Advanced Electrophoresis Solutions (AES) Ltd.**

AES is a company dedicated to providing whole column imaging detection capillary electrophoresis systems, reagents and consumables. AES' total solution encompasses a variety of CEInfinite products and support, which enables clients to focus on protein separation, quantification, and characterization.



Advanced Electrophoresis Solutions Ltd.

1600 Industrial Road  
Cambridge, Ontario N3H 4W5 Canada

Phone: 1 519 653 6888

Fax: 1 519 804 4266

Email: [info@aeslifesciences.com](mailto:info@aeslifesciences.com)

Company Website: [www.aeslifesciences.com](http://www.aeslifesciences.com)

Product Website: [www.ceinfinite.com](http://www.ceinfinite.com)