

# C-Chip

CE

The World's First Genuine  
Disposable Hemocytometer



Precise volume control  
Easy to use &

**NanoEnTek**

e-mail

[sales@nanoentek.com](mailto:sales@nanoentek.com)

website

[www.nanoentek.com](http://www.nanoentek.com)

**NanoEnTek, Inc.**

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,  
Gyeonggi-do, 18531, Korea

Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999

**NanoEnTek America, Inc.**

240 Bear Hill Road, Suite 101, Waltham,  
MA 02451, USA

Tel : +1-781-472-2558 / Fax : +1-781-790-5649

**EC REP**

MT Promedt Consulting GmbH  
Altenhofstrasse 80,  
66386 St. Ingbert, Germany

# C-Chip is



a Disposable Plastic Hemocytometer

- Does not need coverslips
- Eliminates the need of washing and reuse
- Reduces exposures to the infectious samples and contamination
- Guaranteed repeatability and reliability
- Two test
- Quartz grade optical plastic
- Sturdy and strong



## Applications

Blood analysis (Hematology) | Blood cell counting  
 Cell Culture | Cell concentration measurement / Cell viability test  
 Microbiology | Bacteria & Fungal spore counting  
 IVF, IUI | Sperm counting



## Models

No.	Grid	Depth
DHC-N01	Neubauer Improved	100 $\mu\text{m}$
DHC-B01	Burker	100 $\mu\text{m}$
DHC-B02	Burker Turk	100 $\mu\text{m}$
DHC-F01	Fuchs Rosenthal	200 $\mu\text{m}$
DHC-M01	Malassez	200 $\mu\text{m}$

## Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm  
 Chamber depth : 100, 200  $\mu\text{m}$   
 Chamber volume : 10  $\mu\text{L}$ (DHC-N01, B01, B02)  
 20  $\mu\text{L}$ (DHC-F01, M01)

**NanoEnTek**

e-mail  
website

sales@nanointek.com  
www.nanointek.com

**NanoEnTek, Inc.**

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,  
Gyeonggi-do, 18531, Korea  
Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999

**NanoEnTek America, Inc.**

240 Bear Hill Road, Suite 101, Waltham,  
MA 02451, USA  
Tel : +1-781-472-2558 / Fax : +1-781-790-5649

**EC REP**

MT Promedt Consulting GmbH  
Altenhofstrasse 80,  
66386 St. Ingbert, Germany

Your local distributor:



THE INNOVATOR FOR SCIENCE

[WWW.ISOGEN-LIFESCIENCE.COM](http://WWW.ISOGEN-LIFESCIENCE.COM)

**Benelux**

Isogen Life Science B.V.  
Veldzigt 2a  
3454 PW De Meern  
The Netherlands  
tel. +31 30 688 0771

**Spain and Portugal**

Isogen Life Science B.V.  
Veldzigt 2a  
3454 PW De Meern  
The Netherlands  
tel. +34 936 365 178

For service questions, please contact  
[service@isogen-lifescience.com](mailto:service@isogen-lifescience.com)

For all other questions, please contact  
[support@isogen-lifescience.com](mailto:support@isogen-lifescience.com)