



## Centrifuge Tubes

STARLAB's new range of centrifuge tubes are designed for the centrifugation of cells, but can be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of  $10^{-6}$  fulfils the highest demands on sterility, meeting EN 556-1 standard (Sterilisation of Medical Devices).

Maximum centrifugation rate: 12,000 x g.

### Product features

- ▶ Clear, sterile polypropylene tubes with green HDPE screw caps
- ▶ White writing surface and clearly printed graduations
- ▶ Longer length screw caps with sealing ring prevents any leakage
- ▶ Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile (SAL  $10^{-6}$ ) material
- ▶ Certified RNase, DNase, DNA and Pyrogen free
- ▶ Racked tubes supplied in eco-friendly recyclable cardboard racks



ORDERING INFORMATION			
Description	Pack Size	Cat. No.	List £
<b>15 ml Tubes</b>			
15 ml Centrifuge Tubes, Conical, Loose (Sterile)	20 x 25	E1415-0100	£105.15
15 ml Centrifuge Tubes, Conical, Racked (Sterile)	20 x 25	E1415-0700	£111.73
<b>50 ml Tubes</b>			
50 ml Centrifuge Tubes, Conical, Loose (Sterile)	20 x 25	E1450-0100	£105.15
50 ml Centrifuge Tubes, Conical, Racked (Sterile)	20 x 25	E1450-0700	£111.73



Passionate for Science.



www.starlab.co.uk

Your local distributor:



THE INNOVATOR FOR SCIENCE

[WWW.ISOGEN-LIFESCIENCE.COM](http://WWW.ISOGEN-LIFESCIENCE.COM)

Isogen Life Science B.V.  
Veldzigt 2a  
3454 PW De Meern  
The Netherlands  
tel. +31 30 688 0771

For service questions, please contact  
[service@isogen-lifescience.com](mailto:service@isogen-lifescience.com)

For all other questions, please contact  
[support@isogen-lifescience.com](mailto:support@isogen-lifescience.com)