



STARfocus


Serological Pipettes

Discover the possibilities!

STARLAB's new extended range of Serological Pipettes cover seven different volumes and, for the first time, includes a 100 ml pipette. The pipettes are constructed from polystyrene for sterile serological and tissue culture applications. Colour-coding on both the pipettes and the packaging provides easy identification. The black printed, easy to read ascending and descending graduations simplify fluid measuring and dispensing, and all pipettes (excl. 100 ml) feature negative graduations for extra working volume. Perfect for use with all popular pipette controllers, including ErgoOne® FAST.

Product features

- ▶ Accuracy $\pm 2\%$
- ▶ Colour-coded for easy identification
- ▶ Bidirectional graduations provide added applicability
- ▶ Negative graduations (excl. 100 ml) allows additional working volume
- ▶ Supplied with a polyolefin filter plug
- ▶ Individually wrapped in traceable paper-plastic peel-back packs showing lot number and expiry date
- ▶ Specially designed carton with perforated tear-strip doubles as efficient, bench top dispenser
- ▶ Gamma sterilized to SAL 10^{-6} (ISO11137)
- ▶ Certified RNase, DNase and pyrogen free
- ▶ Non-cytotoxic and non-haemolytic latex-free USP Class VI polystyrene
- ▶ ISO 9001 and ISO 13485 manufactured











Convenient, dispenser packaging



Serological Pipettes - Order information




Serological Pipettes				
	Description	Pack Size	Cat. No.	Price
	1.0 ml Serological Pipette, Individually Wrapped (Sterile)	500	E4860-0001	
	2.0 ml Serological Pipette, Individually Wrapped (Sterile)	500	E4860-0002	
	5.0 ml Serological Pipette, Individually Wrapped (Sterile)	200	E4860-0005	
	10 ml Serological Pipette, Individually Wrapped (Sterile)	200	E4860-0010	
	25 ml Serological Pipette, Individually Wrapped (Sterile)	150	E4860-0025	
	50 ml Serological Pipette, Individually Wrapped (Sterile)	100	E4860-0050	
	100 ml Serological Pipette, Individually Wrapped (Sterile)	50	E4860-1100	



Specifications				
Pipette Volume	Colour Coding	Increments	Extra Graduations	Total Volume
1.0 ml	Yellow	0.01 ml	0.3 ml	1.3 ml
2.0 ml	Green	0.01 ml	0.3 ml	2.3 ml
5.0 ml	Blue	0.1 ml	3.0 ml	8.0 ml
10 ml	Orange	0.1 ml	3.0 ml	13 ml
25 ml	Red	0.2 ml	8.0 ml	33 ml
50 ml	Purple	0.5 ml	10 ml	60 ml
100 ml	Pink	1.0 ml	n/a	100 ml



Lot number and expiry date



STARLAB distributor

STARLAB International GmbH

Neuer Höltigbaum 38 • 22143 Hamburg

Tel.: +49 (0)40 675 99 39 0 • Fax: +49 (0)40 675 99 39 20 • info@starlab.de

Passionate for Science.



www.starlabgroup.com

Your local distributor:



THE INNOVATOR FOR SCIENCE

WWW.ISOGEN-LIFESCIENCE.COM

Isogen Life Science B.V.
Veldzigt 2a
3454 PW De Meern
The Netherlands
tel. +31 30 688 0771

For service questions, please contact
service@isogen-lifescience.com

For all other questions, please contact
support@isogen-lifescience.com