

## iFuge D08 (Doctor Centrifuge)

Speed: Max 4500 RPM (1985 g) | Max Capacity: 4 x 50 ml / 8 x 15 ml tubes



Maintenance free  
long life BLDC motor



2 programmable preset  
mode (P1 and P2)



Imbalance detection  
with auto-cut off



Timer Setting 1 to 99 min.  
& infinite mode



Safe low  
voltage device



Lid Lock safety,  
lid opens automatically  
on run completion



### ROTOR OPTIONS



4 x 50ml | 2200 g



8 x 15ml | 2200 g

For more details visit us at [www.neuation.com](http://www.neuation.com)